


Anti-Dnal1 Rabbit Polyclonal Antibody For Wb And Elisa - Q4Ldg9

Item Number: CM0014565



Introduction

High-quality rabbit polyclonal antibody targeting DNAL1, suitable for WB and ELISA applications. Recognizes DNAL1 in human, mouse, and rat samples. Ideal for ciliary function and primary ciliary dyskinesia research.

[Learn More](#)

Parameter	Value
Product Name	DNAL1 Rabbit pAb
Remarks/Alias	LC1, CILD16, C14orf168, DNAL1
Species	Human
GeneID (Human)	83544
GeneID	83544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-190 of human DNAL1 (NP_113615.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q4LDG9
Protein Weight	22kDa
Shipping	Ice bag