


Anti-Dnal1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q4Ldg9

Item Number: CM0024708



Introduction

Rabbit polyclonal antibody targeting DNAL1 (dynein axonemal light chain 1). Validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for ciliary function and PCD studies.

[Learn More](#)

Parameter	Value
Product Name	DNAL1 Rabbit pAb
Remarks/Alias	LC1; CILD16; C14orf168; DNAL1
Species	Human
Gene ID (Human)	83544
Gene ID	83544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-190 of human DNAL1 (NP_113615.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q4LDG9
Protein Weight	22kDa
Shipping	Ice bag