

Anti-Dnase1 Rabbit Polyclonal Antibody For Wb, Elisa - P24855

Item Number: CM0027297



Introduction

High-quality rabbit polyclonal antibody against human DNASE1 (DNase I). Validated for Western blot and ELISA; reacts with human, mouse, and rat samples. Suitable for apoptosis and NET research.

[Learn More](#)

Parameter	Value
Product Name	DNASE1 Rabbit pAb
Remarks/Alias	DNL1; DRNI; DNASE1
Species	Human
GeneID (Human)	1773
GeneID	1773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-282 of human DNASE1 (NP_005214.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24855
Protein Weight	31kDa
Shipping	Ice bag