

Anti-Dnase1L1 Polyclonal Antibody For Wb, If/Icc, Elisa - P49184

Item Number: CM0014486

Introduction

Rabbit polyclonal antibody against human DNASE1L1 (DNase X, DNase I-like 1). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for nuclease and IVD research.

[Learn More](#)

Parameter	Value
Product Name	DNASE1L1 Rabbit pAb
Remarks/Alias	XIB; G4.8; DNL1L; DNASEX; DNAS1L1; DNASE1L1
Species	Human
Gene ID	1774
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human DNASE1L1 (NP_006721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P49184
Protein Weight	34kDa
Shipping	Ice bag