

Anti-Dner Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nft8

Item Number: CM0022295

Introduction

High-quality rabbit polyclonal antibody targeting human DNER. Validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for Notch signaling and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	DNER Rabbit pAb
Remarks/Alias	bet; UNQ26; DNER
Species	Human
GeneID (Human)	92737
GeneID	92737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human DNER (NP_620711.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFT8
Protein Weight	78kDa
Shipping	Ice bag