


Anti-Dnm1 Rabbit Polyclonal Antibody For Wb And Elisa - Q05193

Item Number: CM0023762



Introduction

Rabbit polyclonal antibody targeting human DNM1 (Dynamin-1); validated for Western blot and ELISA; cross-reacts with mouse and rat; ideal for endocytosis and synaptic vesicle recycling research.

[Learn More](#)

Parameter	Value
Product Name	DNM1 Rabbit pAb
Remarks/Alias	DNM, DEE31, EIEE31, DNM1
Species	Human
GeneID (Human)	1759
GeneID	1759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-490 of human DNM1 (NP_004399.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05193
Protein Weight	97kDa
Shipping	Ice bag