

Anti-Dnm2 Polyclonal Antibody For If/Icc, Elisa - P50570

Item Number: CM0015324



Introduction

Rabbit polyclonal antibody targeting human DNM2, suitable for IF/ICC and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying dynamin-2 in endocytosis, actin dynamics, and Charcot-Marie-Tooth disease.

[Learn More](#)

Parameter	Value
Product Name	DNM2 Rabbit pAb
Remarks/Alias	DYN2; CMT2M; DYNII; LCC55; CMTDI1; CMTDIB; DI-CMTB; DNM2
Species	Human
GeneID (human)	1785
GeneID	1785
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 607-866 of human DNM2 (NP_004936.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50570
Protein Weight	98kDa
Shipping	Ice bag