

Anti-Dnm2 Polyclonal Antibody For Wb, If/Icc, Elisa - P50570

Item Number: CM0025956



Introduction

DNM2 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat, targeting the dynamin-2 GTPase involved in endocytosis and cytoskeletal remodeling.

[Learn More](#)

Parameter	Specification
Product Name	DNM2 Rabbit pAb
Remarks/Alias	DYN2; CMT2M; DYNII; LCCS5; CMTDI1; CMTDIB; DI-CMTB; DNM2
Species	Human
GeneID (human)	1785
GeneID	1785
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P50570
Protein Weight	98kDa
Shipping	Ice bag