

Anti-Dnmt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P26358

Item Number: CM0002704



Introduction

DNMT1 Rabbit Monoclonal Antibody validated for WB, IHC-P, IP, ELISA with cross-reactivity to Human, Mouse, Rat. Essential for studying DNA methylation, epigenetic inheritance, and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	DNMT1 Rabbit mAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsa1; DNMT1
Species	Human
GeneID (Human)	1786
GeneID	1786
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DNMT1 (NP_001370.1).
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag