

Anti-Dnmt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P26358

Item Number: CM0016632



Introduction

DNMT1 Rabbit pAb (Human, Mouse, Rat) validated for WB, IHC-P, IF/ICC, IP, ELISA. Targets the maintenance DNA methyltransferase DNMT1 for epigenetic research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	DNMT1 Rabbit pAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsal; DNMT1
Species	Human
Gene ID (Human)	1786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-124 of DNMT1 (NP_001370.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag