


Anti-Dnmt3A Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y6K1

Item Number: CM0002681



Introduction

High-quality rabbit monoclonal antibody against DNMT3A, suitable for WB and ELISA. DNMT3A is a de novo DNA methyltransferase essential for genome-wide methylation patterning during development. Validated for human reactivity, aids epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DNMT3A Rabbit mAb
Remarks/Alias	TBR5; HESJAS; DNMT3A2; M.HsaIIIA; 3A
Species	Human
GeneID (Human)	1788
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human DNMT3A (Q9Y6K1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6K1
Protein Weight	102kDa
Shipping	Ice bag