

Anti-Dnmt3B Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubc3

Item Number: CM0018376



Introduction

Rabbit polyclonal antibody against human DNMT3B, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Essential for de novo DNA methylation during development and gene silencing research.

[Learn More](#)

Parameter	Value
Product Name	DNMT3B Rabbit pAb
Remarks/Alias	ICF; ICF1; FSHD4; M.HsaIIIB; DNMT3B
Species	Human
Gene ID (Human)	1789
Gene ID	1789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human DNMT3B (NP_008823.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBC3
Protein Weight	96kDa
Shipping	Ice bag