

Anti-Dnmt3B Rabbit Monoclonal Antibody For Wb, Chip, Elisa - Q9Ubc3

Item Number: CM0008016



Introduction

Recombinant rabbit monoclonal antibody against human DNMT3B (SWISS Q9UBC3) for Western blot, ChIP, and ELISA applications. DNMT3B is a de novo DNA methyltransferase crucial for establishing methylation patterns during development and silencing repetitive elements.

[Learn More](#)

Parameter	Value
Product Name	DNMT3B Rabbit mAb
Remarks/Alias	ICF; ICF1; FSHD4; M.HsaIIIB
Species	Human
Gene ID (Human)	1789
Gene ID	1789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-360 of human DNMT3B (NP_008823.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ChIP, ELISA
Cross Reactivity	Human
SWISS	Q9UBC3
Protein Weight	96kDa
Shipping	Ice bag