

Anti-Dnmt3B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubc3

Item Number: CM0026465



Introduction

High-quality DNMT3B Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects DNMT3B in human, mouse, rat tissues. Crucial for epigenetic studies of DNA methylation during development.

[Learn More](#)

Parameter	Value
Product Name	DNMT3B Rabbit pAb
Remarks/Alias	ICF; ICF1; FSHD4; M.HsaIIIB; DNMT3B
Species	Human
Gene ID (Human)	1789
Gene ID	1789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 345-475 of human DNMT3B (NP_008823.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBC3
Protein Weight	96kDa
Shipping	Ice bag