


Anti-Dnttip2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Qje6

Item Number: CM0023694



Introduction

High-quality rabbit polyclonal antibody targeting human DNTTIP2 (ERBP/TdIF2). Verified for WB, IF/ICC, ELISA; cross-reacts with mouse. Ideal for transcription regulation and ribosome biogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	DNTTIP2 Rabbit pAb
Remarks/Alias	ERBP; FCF2; TdIF2; HSU15552; LPTS-RP2; DNTTIP2
Species	Human
Gene ID (Human)	30836
Gene ID	30836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human DNTTIP2 (NP_055412.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5QJE6
Protein Weight	84kDa
Shipping	Ice bag