


Anti-Dock2 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q92608

Item Number: CM0027614



Introduction

Anti-DOCK2 rabbit polyclonal antibody validated for WB, IF/ICC, IP, ELISA applications. Recognizes human and mouse DOCK2 (Q92608), a guanine nucleotide exchange factor (GEF) that activates RAC1/RAC2 and regulates lymphocyte migration.

[Learn More](#)

Parameter	Value
Product Name	DOCK2 Rabbit pAb
Remarks/Alias	IMD40; DOCK2
Species	Human
Gene ID (Human)	1794
Gene ID	1794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1551-1830 of human DOCK2 (NP_004937.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92608
Protein Weight	212kDa
Shipping	Ice bag