

Anti-Dock3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Izd9

Item Number: CM0023776

Introduction

Rabbit polyclonal antibody targeting human DOCK3 (PBP/MOCA). Suitable for Western blot and ELISA. Explore Alzheimer's disease mechanisms and cytoskeletal signaling with this validated reagent.

[Learn More](#)

Parameter	Value
Product Name	DOCK3 Rabbit pAb
Remarks/Alias	PBP; MOCA; NEDIDHA; DOCK3
Species	Human
Gene ID	1795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1348-1673 of human DOCK3 (NP_004938.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IZD9
Protein Weight	233KDa
Shipping	Ice bag