


Anti-Dock7 Polyclonal Antibody For Wb, Elisa - Q96N67

Item Number: CM0019599



Introduction

Rabbit polyclonal antibody against human DOCK7 (Q96N67). Suitable for WB and ELISA; reacts with human, mouse, and rat. DOCK7 is a Rac1/3 GEF involved in neuronal polarization and actin cytoskeleton regulation.

[Learn More](#)

Parameter	Value
Product Name	DOCK7 Rabbit pAb
Remarks/Alias	ZIR2; DEE23; EIEE23; DOCK7
Species	Human
GeneID (Human)	85440
GeneID	85440
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human DOCK7 (NP_001354490.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96N67
Protein Weight	243kDa
Shipping	Ice bag