

Anti-Dog-1/Tmem16A Rabbit Polyclonal Antibody For Wb, Elisa - Q5Xxa6

Item Number: CM0020442

Introduction

Rabbit polyclonal antibody against human DOG-1/TMEM16A (ANO1), validated for WB and ELISA with Human and Mouse cross-reactivity. Suitable for studying calcium-activated chloride channels, gastrointestinal motility, and airway secretion.

[Learn More](#)

Parameter	Value
Product Name	DOG-1/TMEM16A Rabbit pAb
Remarks/Alias	DOG1; INDMS; TAOS2; ORAOV2; TMEM16A; ANO1
Species	Human
GeneID (Human)	55107
GeneID	55107
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-210 of human ANO1 (NP_060513.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5XXA6
Protein Weight	114kDa
Shipping	Ice bag