


Anti-Dok1 Monoclonal Antibody For Wb, If/Icc, Elisa - Q99704

Item Number: CM0002319



Introduction

Rabbit monoclonal antibody targeting human DOK1 (docking protein 1). Validated for WB, IF/ICC, and ELISA applications with cross-reactivity to human and mouse samples. Useful for studying insulin signaling and integrin activation pathways.

[Learn More](#)

Parameter	Value
Product Name	DOK1 Rabbit mAb
Remarks/Alias	pp62; P62DOK; DOK1
Species	Human
GeneID (Human)	1796
GeneID	1796
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 382-481 of human DOK1 (Q99704)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99704
Protein Weight	52kDa
Shipping	Ice bag