


Anti-Dok1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99704

Item Number: CM0027825



Introduction

Rabbit polyclonal antibody against DOK1 (Docking protein 1, p62(dok), pp62). Validated in human, mouse, rat for WB, IHC-P, IF/ICC, ELISA. Targets insulin signaling and integrin activation pathways.

[Learn More](#)

Parameter	Value
Product Name	DOK1 Rabbit pAb
Remarks/Alias	pp62; P62DOK; DOK1
Species	Human
GeneID (Human)	1796
GeneID	1796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DOK1 (NP_001372.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99704
Protein Weight	52kDa
Shipping	Ice bag