

Anti-Dok4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Tew6

Item Number: CM0016560

Introduction

Rabbit polyclonal antibody targeting human DOK4 (IRS5). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for studying RET signaling and neuronal differentiation.

[Learn More](#)

Parameter	Specification
Product Name	DOK4 Rabbit pAb
Remarks/Alias	IRS5; IRS-5; DOK4
Species	Human
Gene ID	55715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human DOK4 (NP_060580.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEW6
Protein Weight	37kDa
Shipping	Ice bag