


Anti-Dok5 Polyclonal Antibody For Wb, Elisa - Q9P104

Item Number: CM0018604



Introduction

Rabbit polyclonal antibody against human DOK5 (IRS-6), validated for Western blotting and ELISA. Cross-reacts with mouse and rat. Suitable for studies of RET signaling, neuronal differentiation, and MAP kinase pathway.

[Learn More](#)

Parameter	Value
Product Name	DOK5 Rabbit pAb
Remarks/Alias	IRS6; IRS-6; C20orf180; DOK5
Species	Human
GeneID (human)	55816
GeneID	55816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 237-306 of human DOK5 (NP_060901.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9P104
Protein Weight	35kDa
Shipping	Ice bag