


Anti-Dopamine Beta-Hydroxylase (Dbh) Rabbit Polyclonal Antibody For Wb, Elisa - P09172

Item Number: CM0020404



Introduction

Rabbit polyclonal antibody targeting human DBH, suitable for WB and ELISA, with cross-reactivity to mouse and rat. Enables study of dopamine hydroxylation to noradrenaline.

[Learn More](#)

Parameter	Specification
Product Name	DBH Rabbit pAb
Remarks/Alias	DBM, ORTHYP1, DBH
Species	Human
Gene ID (Human)	1621
Gene ID	1621
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-200 of human DBH (NP_958852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09172
Protein Weight	69kDa
Shipping	Ice bag