


Anti-Dopamine Beta Hydroxylase Rabbit Monoclonal Antibody For Wb And Elisa - P09172

Item Number: CM0007298



Introduction

High-quality rabbit monoclonal antibody targeting human Dopamine beta Hydroxylase (DBH), suitable for WB and ELISA. Essential tool for studying catecholamine biosynthesis, neurotransmitter regulation, and related disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Dopamine beta Hydroxylase Rabbit mAb
Remarks/Alias	DBM; ORTHYP1; Dopamine beta Hydroxylase
Species	Human
GeneID (Human)	1621
GeneID	1621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human Dopamine beta Hydroxylase (NP_000778.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09172
Protein Weight	69kDa
Shipping	Ice bag