

Anti-Dopamine Receptor D3 (Drd3) Polyclonal Antibody For Wb, Elisa - P35462

Item Number: CM0017809



Introduction

High-quality rabbit polyclonal antibody targeting human Dopamine Receptor D3 (DRD3) for WB and ELISA. Suitable for detection in mouse and rat samples. For research on dopamine signaling, neuropsychiatric disorders, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	Dopamine Receptor D3 (DRD3) Rabbit pAb
Remarks/Alias	D3DR; ETM1; FET1; Dopamine Receptor D3 (DRD3)
Species	Human
GeneID	1814
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human DRD3 (NP_000787.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35462
Protein Weight	44kDa
Shipping	Ice bag