

Anti-Dopamine Receptor D3 (Drd3) Rabbit Monoclonal Antibody For Wb, Elisa - P35462

Item Number: CM0008958



Introduction

Rabbit monoclonal antibody against human Dopamine Receptor D3 (DRD3), for WB and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for neuroscience research on dopamine signaling and neuropsychiatric disorders.

[Learn More](#)

Parameter	Value
Product Name	Dopamine Receptor D3 (DRD3) Rabbit mAb
Remarks/Alias	D3DR; ETM1; FET1; Dopamine Receptor D3 (DRD3)
Species	Human
GeneID (Human)	1814
GeneID	1814
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Dopamine Receptor D3 (DRD3) (P35462).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35462
Protein Weight	44kDa
Shipping	Ice bag