

Anti-Dot1L/Kmt4 Polyclonal Antibody For Wb, Elisa - Q8Tek3

Item Number: CM0017729



Introduction

Rabbit polyclonal antibody targeting human Dot1L/KMT4 (H3K79 methyltransferase), validated in WB and ELISA, cross-reacts with rat and mouse. Ideal for studying chromatin regulation and gene silencing.

[Learn More](#)

Parameter	Value
Product Name	Dot1L/KMT4 Rabbit pAb
Remarks/Alias	DOT1; KMT4; Dot1L/KMT4
Species	Human
Gene ID (Human)	84444
Gene ID	84444
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Dot1L/KMT4 (NP_115871.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat, Mouse
SWISS	Q8TEK3
Protein Weight	165kDa
Shipping	Ice bag