

Rabbit polyclonal antibody targeting human Doublecortin (DCX), validated for Western blot, immunofluorescence/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for neuronal migration and cortical development studies.

Parameter	Value
Product Name	Doublecortin Rabbit pAb
Remarks/Alias	DC; DBCN; LISX; SCLH; XLIS; Doublecortin
Species	Human
Gene ID (Human)	1641
Gene ID	1641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-365 of human Doublecortin (NP_835365.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43602
Protein Weight	41kDa
Shipping	Ice bag