

Anti-Doublecortin Rabbit Polyclonal Antibody  
validated for WB, IF/ICC, IF-P, ELISA; reacts with  
human, mouse, rat. Manufactured using  
recombinant immunogen; ideal for neuronal  
migration and cortical lamination studies.

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Doublecortin Rabbit pAb                                                                                                   |
| Remarks/Alias    | DC; DBCN; LISX; SCLH; XLIS; Doublecortin                                                                                  |
| Species          | Human                                                                                                                     |
| GenelD (Human)   | 1641                                                                                                                      |
| GenelD           | 1641                                                                                                                      |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 265-365 of human Doublecortin (NP_835365.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | O43602                                                                                                                    |
| Protein Weight   | 41kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |