


Anti-Doublecortin Rabbit Monoclonal Antibody For Wb, Elisa - 043602

Item Number: CM0006384



Introduction

Doublecortin Rabbit monoclonal antibody validated for Western blot and ELISA applications, detecting human, mouse, and rat Doublecortin (41 kDa). Suitable for neurodevelopment and neuronal migration studies.

[Learn More](#)

Parameter	Value
Product Name	Doublecortin Rabbit mAb
Remarks/Alias	DC; DBCN; LISX; SCLH; XLIS; Doublecortin
Species	Human
Gene ID (Human)	1641
Gene ID	1641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-365 of human Doublecortin. (NP_835365.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43602
Protein Weight	41kDa
Shipping	Ice bag