

Doublecortin Rabbit Monoclonal Antibody
validated for WB and ELISA with cross-reactivity
to mouse and rat. Ideal for neuronal migration
and developmental biology research.

Parameter	Value
Product Name	Doublecortin Rabbit mAb
Remarks/Alias	DC, DBCN, LISX, SCLH, XLIS, Doublecortin
Species	Human
Gene ID (Human)	1641
Gene ID	1641
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Doublecortin (O43602)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43602
Protein Weight	41kDa
Shipping	Ice bag