


Anti-Dp1/Tfdp1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14186

Item Number: CM0026491



Introduction

Rabbit polyclonal antibody against human DP1/TFDP1 (Transcription factor Dp-1). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Ideal for cell cycle and transcription research.

[Learn More](#)

Parameter	Value
Product Name	DP1/TFDP1 Rabbit pAb
Remarks/Alias	DP1; DILC; Dp-1; DRTF1; DP1/TFDP1
Species	Human
GeneID (Human)	7027
GeneID	7027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human DP1/TFDP1 (NP_009042.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14186
Protein Weight	45kDa
Shipping	Ice bag