


Anti-Dp1/Tfdp1 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q14186

Item Number: CM0013111



Introduction

Rabbit polyclonal antibody targeting human DP1/TFDP1, a transcription factor involved in E2F-dependent cell cycle regulation. Validated for IHC-P and ELISA with cross-reactivity to mouse and rat. ~45 kDa.

[Learn More](#)

Parameter	Value
Product Name	DP1/TFDP1 Rabbit pAb
Remarks/Alias	DP1; DILC; Dp-1; DRTF1; DP1/TFDP1
Species	Human
GeneID (Human)	7027
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human DP1/TFDP1 (NP_009042.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14186
Protein Weight	45kDa
Shipping	Ice bag