


Anti-Dpep1 Rabbit Polyclonal Antibody For Wb, Elisa - P16444

Item Number: CM0025070



Introduction

DPEP1 Rabbit Polyclonal Antibody validated for WB and ELISA, recognizes human and mouse DPEP1 (~46 kDa). DPEP1 hydrolyzes dipeptides and serves as an adhesion receptor for neutrophils. Sourced from rabbit, ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	DPEP1 Rabbit pAb
Remarks/Alias	MDP; RDP; MBD1; DPEP1
Species	Human
GeneID (Human)	1800
GeneID	1800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-385 of human DPEP1 (NP_001121613.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16444
Protein Weight	46kDa
Shipping	Ice bag