

## Anti-Dpep2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H4A9

Item Number: CM0019441

### Introduction

Rabbit polyclonal antibody against human DPEP2, validated for WB and ELISA. Reacts with human, mouse, and rat. DPEP2 is a dipeptidase involved in leukotriene metabolism and inflammatory signaling.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | DPEP2 Rabbit pAb                                                                                                    |
| Remarks/Alias    | MBD2; DPEP2                                                                                                         |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 64174                                                                                                               |
| GeneID           | 64174                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 400-463 of human DPEP2 (NP_071750.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9H4A9                                                                                                              |
| Protein Weight   | 53kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |