

Anti-Dpf2 Rabbit Polyclonal Antibody For Wb And Elisa - Q92785

Item Number: CM0023548



Introduction

High-quality rabbit polyclonal antibody against human DPF2, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for chromatin research and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	DPF2 Rabbit pAb
Remarks/Alias	REQ, CSS7, UBID4, ubi-d4, SMARCG2, DPF2
Species	Human
Gene ID (Human)	5977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human DPF2 (NP_006259.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92785
Protein Weight	44kDa
Shipping	Ice bag