


Anti-Dpf3 Rabbit Polyclonal Antibody For Wb, Elisa - Q92784

Item Number: CM0021590



Introduction

High-quality rabbit polyclonal antibody targeting human DPF3 (also known as BAF45C). Validated for WB and ELISA, with cross reactivity to mouse and rat. Suitable for chromatin remodeling and neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DPF3 Rabbit pAb
Remarks/Alias	CERD4, BAF45C, SMARCG3, DPF3
Species	Human
GeneID (human)	8110
GeneID	8110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human DPF3 (NP_036206.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92784
Protein Weight	43kDa
Shipping	Ice bag