

Anti-Dph2 Polyclonal Antibody For Wb, Elisa - Q9Bqc3

Item Number: CM0022095

Introduction

Rabbit polyclonal antibody targeting human DPH2, validated for WB and ELISA. Detects mouse DPH2 cross-reactively. Suitable for diphthamide biosynthesis and elongation factor research.

[Learn More](#)

Parameter	Value
Product Name	DPH2 Rabbit pAb
Remarks/Alias	DPH2L2; DEDSSH2; DPH2
Species	Human
GeneID (Human)	1802
GeneID	1802
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 120-330 of human DPH2 (NP_001375.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BQC3
Protein Weight	52kDa
Shipping	Ice bag