

Anti-Dpm1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60762

Item Number: CM0026242

Introduction

High-quality rabbit polyclonal antibody to DPM1, the catalytic subunit of DPM synthase, validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for studying N-glycosylation and GPI anchor synthesis.

[Learn More](#)

Parameter	Value
Product Name	DPM1 Rabbit pAb
Remarks/Alias	MPDS; CDGIE; DPM1
Species	Human
GeneID (human)	8813
GeneID	8813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DPM1 (NP_003850.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60762
Protein Weight	30kDa
Shipping	Ice bag