

Anti-Dpp3 Polyclonal Antibody For Wb, Elisa - Q9Ny33

Item Number: CM0017247



Introduction

Recombinant protein-immunized rabbit polyclonal antibody against human DPP3, suitable for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying dipeptidyl peptidase III functions in peptide degradation.

[Learn More](#)

Parameter	Value
Product Name	DPP3 Rabbit pAb
Remarks/Alias	DPPIII; DPP3
Species	Human
GeneID (human)	10072
GeneID	10072
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DPP3 (NP_569710.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY33
Protein Weight	83kDa
Shipping	Ice bag