

Anti-Dpp3 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Ny33

Item Number: CM0006030



Introduction

Rabbit monoclonal antibody against human DPP3, validated for Western blot (WB) and ELISA. Detects DPP3 in human, mouse, and rat. Ideal for research on bioactive peptide degradation.

[Learn More](#)

Parameter	Value
Product Name	DPP3 Rabbit mAb
Remarks/Alias	DPPIII; DPP3
Species	Human
Gene ID (Human)	10072
Gene ID	10072
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 638-737 of human DPP3 (Q9NY33).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY33
Protein Weight	83kDa
Shipping	Ice bag