

Anti-Dpp8 Polyclonal Antibody For Wb And Elisa - Q6V1X1

Item Number: CM0026402



Introduction

Rabbit polyclonal antibody targeting human DPP8, raised against the N-terminal region (aa10-230). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for pyroptosis and inflammasome studies.

[Learn More](#)

Parameter	Value
Product Name	DPP8 Rabbit pAb
Remarks/Alias	DP8; DPRP1; DPRP-1; MST097; MSTP097; MSTP135; MSTP141; DPP8
Species	Human
Gene ID (Human)	54878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-230 of human DPP8 (NP_932065.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6V1X1
Protein Weight	103kDa
Shipping	Ice bag