

Anti-Dppa2 Polyclonal Antibody For Wb, Elisa - Q7Z7J5

Item Number: CM0015063

Introduction

Rabbit polyclonal antibody against DPPA2, validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for pluripotency and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	DPPA2 Rabbit pAb
Remarks/Alias	CT100; PESCRG1; ECAT15-2; DPPA2
Species	Human
GeneID (Human)	151871
GeneID	151871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DPPA2 (NP_620170.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z7J5
Protein Weight	34kDa
Shipping	Ice bag