


Anti-Dppa3/Stella Polyclonal Antibody For Wb And Elisa - Q8Qzy3

Item Number: CM0019760



Introduction

Rabbit polyclonal antibody targeting DPPA3/Stella for WB and ELISA. Recognizes mouse and rat Stella; key regulator of epigenetic reprogramming and maternal DNA methylation protection in early embryos. Ideal for developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Stella Rabbit pAb
Remarks/Alias	PCG7, PGC7, Stella, 2410075G02Rik
Species	Mouse
Gene ID (Human)	73708
Immunogen	Recombinant fusion protein containing human Stella (aa 1-150)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8QZY3
Protein Weight	18kDa
Shipping	Ice bag