


Anti-Dpys Polyclonal Antibody For Wb, Elisa - Q14117

Item Number: CM0025780



Introduction

Rabbit polyclonal antibody against human DPYS (dihydropyrimidinase), validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for pyrimidine metabolism and related disease research.

[Learn More](#)

Parameter	Value
Product Name	DPYS Rabbit pAb
Remarks/Alias	DHP; DHPase; DPYS
Species	Human
GeneID (human)	1807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-519 of human DPYS (NP_001376.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14117
Protein Weight	57kDa
Shipping	Ice bag