


Anti-Dpysl3 Rabbit Polyclonal Antibody For Wb, Elisa - Q14195

Item Number: CM0016730



Introduction

High-quality rabbit polyclonal antibody against human DPYSL3 (CRMP-4). Suitable for WB and ELISA, cross-reacts with mouse and rat. Ideal for axon guidance and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	DPYSL3 Rabbit pAb
Remarks/Alias	DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1; DPYSL3
Species	Human
GeneID	1809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-520 of human DPYSL3 (NP_001378.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14195
Protein Weight	62kDa
Shipping	Ice bag