

## Anti-Dr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q01658

Item Number: CM0028458



### Introduction

High-quality anti-DR1 rabbit polyclonal antibody validated for WB and ELISA. Detects DR1 (NC2/NC2-beta) across human, mouse, and rat samples. Ideal for transcription regulation research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | DR1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | NC2; NC2B; NCB2; NC2-BETA; DR1                                                                                  |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 1810                                                                                                            |
| GeneID           | 1810                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human DR1 (NP_001929.1). |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | Q01658                                                                                                          |
| Protein Weight   | 19kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |