


Anti-Dr4 Rabbit Monoclonal Antibody For Fc, Elisa - 000220

Item Number: CM0007533



Introduction

Rabbit monoclonal antibody targeting human DR4 (TNFRSF10A/TRAIL-R1), validated for flow cytometry (FC) and ELISA. Recombinant immunogen (aa 124-232) ensures specificity. Ideal for apoptosis and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	DR4 Rabbit mAb
Remarks/Alias	TNFRSF10A; APO2; CD261; DR4; TRAILR-1; TRAILR1; TNF receptor superfamily member 10a
Species	Human
GeneID (Human)	8797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-232 of human DR4 (NP_003835.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	O00220
Protein Weight	50kDa
Shipping	Ice bag