


Anti-Dr5/Trailr2/Tnfrsf10B Rabbit Polyclonal Antibody For Wb And Elisa - O14763

Item Number: CM0015161



Introduction

Rabbit polyclonal antibody targeting DR5/TRAILR2/TNFRSF10B, validated for Western blot and ELISA. Specific for human death receptor 5, a key mediator of TRAIL-induced apoptosis. Ideal for apoptosis and cancer research. 48kDa.

[Learn More](#)

Parameter	Value
Product Name	DR5/TRAILR2/TNFRSF10B Rabbit pAb
Remarks/Alias	DR5; CD262; KILLER; TRICK2; TRICKB; ZTNFR9; TRAILR2; TRICK2A; TRICK2B; TRAIL-R2; KILLER/DR5
Species	Human
GeneID (human)	8795
GeneID	8795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-440 of human DR5/TRAILR2/TNFRSF10B (NP_003833.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14763
Protein Weight	48kDa
Shipping	Ice bag